



Actual Size:
148 mm (Width)
x 110 mm (Height)
or
5³/₄" (Width)
x 4¹/₄" (Height)




Ruler is for scaling purposes.


How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

AirTouch™ **N30i**

Nasal cradle mask





ResMed

AirTouch™ N30i


AirFit™ N30i

Nasal mask

AirFit™ P30i

Nasal pillows mask

Frame Fitting Instructions



 Small frame	 Standard frame
0–150mm	>150mm
0–6"	>6"